## **XR75 Seal Check**





High-sensitivity seal inspection ideal for patches, transdermal, and oral film products

Model: KXE7510DGEKE

- Complies with FDA 21 CFR Part 11
- Eliminates packages with product trapped in seal area
- Great for pharmaceutical & nutraceutical applications
- Also includes contaminant and product defect inspection (e.g. missing, broken, size)
- Split conveyor facilitates highest sensitivity
- Operator-friendly product setup
- Packaged product applications

## **SPECS**

<b>Overall Length</b>	800mm
Max Product Height	50mm
Max Product Width	205mm



